


## Abstract

Traditional SPR-based drug discovery efforts in the small-molecule space have revolved around single-use chips for screening and hit-to-lead efforts. Typically, proteins are recombinantly or chemically biotinylated for capture on custom-made or, more recently, pre-coated (strept)avidin chips. These surfaces are unable to be regenerated for reuse with freshly immobilized target protein due to the sub-picomolar affinity of the biotin-avidin interaction and the durability of the capture-protein to harsh chemical exposure. Inspired by recent technological advances, we have established simple regenerable protocols that can (1) accommodate target-proteins with common tags such as His- and Avi-, (2) reach sufficient surface density for small-molecule testing, and (3) maintain a stable-baseline for accurate kinetic measurements. On a Biacore 8k+ system, these protocols make effective use of chips, saving time and reducing waste. Additionally, the simple replenishment process allows for special applications such as testing challenging targets with poor longevity on the surface, or quasi-stable multi-component complexes. Through regeneration, we can efficiently test small covalent ligands, a process that was previously considered low-throughput and expensive. The methods are broadly applicable for both early-stage screening and late-stage characterizations by measuring affinities of candidates with long residence times, leveraging residual occupancy (aka chaser) formats.



## Comparison of key SPR-based capture strategies

	Amine Coupling	Avidin-Biotin	Oligo-avidin	NTA (Ni-affinity)	His-avidin	Antibody
Chip Type	CM5	SA / NA	CAP / RGD	NTA	NTA	CM5 / PAG
Regenerable	X	X	✓	✓	✓	✓
Density	High	High	Low/medium	High	Medium	Medium/low
Stability	+++	+++	+++	+ / ++	++ / +++	++
Oriented	No	Yes	Yes	Yes	Yes	Yes
PTS	Moderate	High	Moderate	Moderate	High	Moderate



## Model System for Small-molecules: Carbonic Anhydrase



## Strategy 1: His-Streptavidin (NTA Chip)





## Strategy 2: Oligo-Streptavidin (XanTec RGD chip)



## Advanced Applications: Chaser, Ternary complex

Highly stable and regenerable surface opens opportunities to measure long off-rates. Sufficient stability observed on His-Streptavidin surface:



## Summary

- 2 highly effective capture strategies have been implemented at WuXi AppTec to perform sophisticated SPR analysis of small molecules using reproducible, reusable chip surfaces
- His-Streptavidin makes highly stable complexes, allowing for direct testing of small molecules
  - Can regenerate to provide high-quality kinetic measurements or apply multiple walk-away cycles of advanced formats such as the chaser
- A functionally similar oligo-streptavidin surface can accomplish the same tasks, though with slightly lower protein densities

## References

Gunnarsson A, Stubbs C, Rawlins P, Taylor-Newman E, Lee W, Geschwindner S, Hytönen V, Holdgate G, Jha R, Dahl G. Regenerable Biosensors for Small-Molecule Kinetic Characterization Using SPR. *SLAS Discovery*. 2021; 26 (5): 730-739

Papalia G, Leavitt S, Bynum M, Katsamba P, Wilton R, Qiu H, Steukers M, Wang S, Bindu L, Phogat S, Gianetti A, Ryan T, Pudlak V, Matusiewicz K, Michelson K, Nowakowski A, Pham-Baginski A, Brooks J, Tieman B, Bruce B, Vaughn M, Baksh M, Cho Y, De Wit M, Smets A, Vandersmissen J, Michiels L, Myszka D. Comparative analysis of 10 small molecules binding to carbonic anhydrase II by different investigators using Biacore technology. *Anal Biochem*. 2006 359 (1): 94-105

Quinn J. Evaluation of Taylor dispersion injections: determining kinetic/affinity interaction constants and diffusion coefficients in label-free biosensing. *Anal Biochem* 2012 Feb 15;421(2):401-10

Zollman D, Shepherd C, Cullis A. Cytiva Application Note, Characterizing PROTAC ternary complex formation using Biacore™ SPR systems

Quinn J, Pitts K, Steffek M, and Mulvihill M, Determination of Affinity and Residence Time of Potent Drug-Target Complexes by Label-free Biosensing. *J Med Chem* 2018 61 (12), 5154-5161