

On-bead DEL's Application and Case Study on Drug Discovery

WuXi Biology

Xiaobing Zhang, Qiaoqiao Zhu, Mengtian Huang, Jason Deng, Yage Liang, Chaoyang Zhang, Meng Guo, Zixuan Zhang, Yihui Xie, Siqi Qian, Chang Liu, Alex Satz, Zhifeng Yu, Zhongyao Ma, Wenji Su, Letian Kuai
WuXi AppTec Headquarters, 288 Fute Zhong Road Waigaoqiao Free Trade Zone, Pudong District, Shanghai, 200131, China.

Introduction

DNA-encoded library (DEL) technology has been recognized as one of the major screening methods with its unique advantages in vast chemical diversity and multiplexed affinity-based screening. Recently, solid-phase On-bead DEL, also known as one bead one compound (OBOC), has been developed to further expand DEL screening from affinity-based screening into biochemical activity screening. Here, we present multiple application case study of OBOC, including biochemical activity screening and traditional affinity screening. In biochemical screening, with high-encapsulation rate microfluidics chip design, we successfully screened a 200K-member OBOC library. Further validation shows that potent inhibitor can be directly identified. In affinity screening, we show that molecular glue can be screened using special focus library along with affinity-based screening methods. Our study demonstrates that OBOC technology could be used for wide-range of drug discovery scenarios.






Figure 1. The schematic diagram of On-bead DEL

Chemistry Reaction & On-Bead DEL Library

On-Bead DEL Reaction Tool Box			
Type	Reaction	Type	Reaction
Amine	Amide formation	Aldehyde	Reductive amination
	Urea		Knoevenagel
	Thiourea	Alcohol	Mitsunobu
	Carbamate	Reduction	Azide reduction
	Sulfonylation	Macrocyclization	Cyclic Peptoid
	Tosylation	Hydrolysis	Methyl/Ethyl ester
Coupling	S _N Ar	Cyclization	Van Leusen for imidazole
	SN2		Click for triazole
	Chan-lam		Oxadiazole formation
	Michael addition		Thiadiazole formation
	Curtius rearrangement	De-protection	De-Fmoc, De-Boc, De-Tfa, De-Trt
Coupling	Guanidination		Heck
	Suzuki	Under development	Isoxazole/isoxazoline
	Sonogashira		Indole
Coupling	Buchwald		Hayashi Miyaura
	Wittig		



Primarily 2-cycle DEL	
<ul style="list-style-type: none">Privileged and target-specific scaffoldsExpected library size: 50-300KPhotocleavable linker: o-nitrobenzene or coumarine-based core, reacted group (-X, -OH, -NH2, -CO)Diverse leaving groups in released compounds	
Drug-like Properties	
<ul style="list-style-type: none">MW ≤500 g/mol, logP ≤5, HBD ≤5, HBA ≤10, PSA ≤140 ÅNo PAINS, No cytotoxic structures	



On-bead DELs Enable Molecular Glue (MG) Discovery




Figure 2. The schematic diagram of CRBN focused DEL and FACS-based glue screening (A) The schematic design of on-bead DEL. The on-bead CRBN focus library is off-the shelf. (B) The overall principle of on-bead DEL MG screening. Following complex formation, protein staining is performed using corresponding fluorescence-conjugated antibodies. Beads exhibiting positive fluorescence signals are subsequently sorted out by FACS.






Figure 3. On-bead CRBN focused library screening identified novel molecular glue (A) The representative FACS sorting figures in the indicated condition during on-bead CRBN focus library screening. (B) The total bead count and recovery rate of spike-in reference compound across different post-selection samples after NGS sequencing. (C) Structures of the five selected compounds. Compound 3 and 4 highlighted in blue were successfully validated. (D) Validation results of hit compounds by HTRF assay. IC50 of compound 3 and 4 are at nM level. (E) Representative validation results of off-DNA compound 4 by SPR. KD value is 70.6 nM. (F) Cellular degradation assay of GSPT1 in KG-1 after 24-hour treatment with the indicated concentrations of CC 90009 or compound 3 or 4.

Activity-Based On-bead DEL Screen for Enzyme Inhibitor



Capability of On-bead DEL on cellular assay



Summary

1. Solid-phase synthesis enables reactions and hot core scaffolds inaccessible to traditional DEL. Build 2-cycle 100K library, expect timeline: 4-6 weeks
2. Affinity based screening capability. Combined with FACS-based sorting, On-bead DEL facilitate molecular glue drug discovery by enabling the direct detection of the ternary complex. Screening Library size <500K, 30x library, expect timeline: 8 weeks for final report (pre-selection, selection, data analysis)
3. High encapsulation chip enables higher throughput of on-Bead DEL screening and cellular on-bead DEL assays. Fast & accurate sorting system with sorting purity higher than 95%. On-Bead DEL specialized NGS & post selection data analysis pipeline. Mpro Assay based on-bead DEL selection shows > 28% validation success rate. High-throughput screening enables 2 day screening and sorting of library size 100 K, 20X library.

