

Evaluations of Targeted Delivery Systems for RNA Therapeutics

Yurun Zhang¹, Dihang Lin¹, Gang Liu¹, Shunchuan Zhang¹, Andy Ng¹, Wenji Su¹

¹Early Discovery Platform, WuXi Biology, WuXi AppTec

WuXi Biology

Introduction

- RNA therapeutics can target or deliver protein-encoding mRNAs with high specificity, an advantage to treat diseases with known genetic mutations. The development of RNA therapies is hindered by a lack of efficient and specific delivery strategy. In recent years, significant progress has been made in conjugate-based and carrier-based delivery systems to improve safe and targeted delivery of oligonucleotides to a variety of tissues and organs.
- To effectively evaluate the capability of novel delivery systems, we have established a set of *in vitro* and *ex vivo* assays to assess the cellular uptake and functional delivery by novel delivery strategies. Here, we demonstrate our capability to assess the binding and internalization of GalNAc conjugates; and to visualize and quantify the cellular uptake, endosomal escape, and functional delivery of lipid nanoparticle (LNP) delivery.
- To overcome the liver tropism of LNPs and expand its therapeutic indications, active-targeted LNPs becomes a strategy for extrahepatic delivery. Here, we generated two types of active targeted LNPs through co-formulation with PDL1-targeting peptide and via surface conjugation with anti-CD22 antibody for specific tissue delivery.

Methods

Flow Cytometric Detection of binding and internalization of GalNAc conjugates


To measure the binding between GalNAc conjugates and ASGPR: Hepatic cells were incubated with GalNAc-Cy5 for 4 hours on ice followed by quantification of fluorescence intensities by flow cytometry. Regarding cellular uptake of GalNAc-Cy5 or GalNAc-conjugated siRNA compounds, hepatic cells were treated with test compounds for 4 hours at 37 °C. Treated cells were then harvested for flow cytometric analysis. To evaluate the uptake efficiency of unlabeled GalNAc conjugates: Cy5-labeled GalNAc (EE) was mixed with GalNAc (L96) conjugates at different concentrations and then incubated with hepatocytes for 4 – 6 hours. Following incubation, cells were collected to measure the endocytosis of GalNAc (EE) -Cy5 via flow cytometry.

Evaluating the cellular uptake, endosomal escape and functional delivery of LNP


HeLa and HeLa-Gal9-GFP cells were treated with MC3-LNP containing Cy5-labeled, GFP-encoding mRNAs across a dose range of 0.1 – 1 µg/mL. Cellular uptake (Cy5) and functional delivery (GFP) were evaluated by measuring the fluorescence intensities. At 4, 6, and 16 hours post treatment, LNP-treated cells were harvested for flow cytometric analysis and confocal microscopy imaging. For live-cell imaging, MC3-LNP treated cells were imaged every 4 hours for 48 hours. To evaluate the endosomal escape of tested LNPs, HeLa-Gal9-GFP cells were imaged live using Yokogawa CQ1 high-content confocal imaging system. Foci formation by Gal9 was counted to quantify the level of endosomal escape followed by LNP treatment *in vitro*.

Results


Binding and internalization of GalNAc conjugates




Visualization and quantification of cellular uptake and delivery of LNP via live-cell imaging, confocal microscopy, and flow cytometry



Active-targeted LNP as a strategy for extrahepatic delivery



Utilizing Galectin-9 foci formation to monitor endosomal escape of LNP



Conclusions

We have established a set of assays to evaluate the uptake and delivery of conjugate-based and carrier-based delivery systems, including: 1) The binding and uptake of GalNAc conjugates into primary hepatocytes; 2) The cellular uptake and delivery of mRNA by LNP; 3) the endosomal escape of LNP in Galectin-9 cell line. In addition, we also demonstrated two strategies to generate active targeted LNPs via 1) co-formulation of PDL1 targeting peptide and 2) conjugation of anti-CD22 antibody onto LNPs for specific organ targeting. In summary, our team has accumulated extensive experiences in assessing *in vitro* delivery of oligonucleotides conjugated with various delivery molecules. Our results also demonstrate our capability to generate active targeted LNPs via post-modification and co-formulation. The platform we have established will further empower potential clients to accelerate discovery and development process of RNA therapeutics.

References

- Paunovska K, Loughrey D, Dahlman J E. Drug delivery systems for RNA therapeutics[J]. Nature Reviews Genetics, 2022, 23(5): 265-280.
- Munson M J, O'Driscoll G, Silva A M, et al. A high-throughput Galectin-9 imaging assay for quantifying nanoparticle uptake, endosomal escape and functional RNA delivery[J]. Communications biology, 2021, 4(1): 211.
- Kornilova, Anna Y , et al. Development of a fluorescence polarization binding assay for asialoglycoprotein receptor. Analytical Biochemistry 425.1(2012):43-46
- Michael T , Lisa H , Alexey R ,et al. Asialoglycoprotein receptor 1 mediates productive uptake of N-acetylgalactosamine-conjugated and unconjugated phosphorothioate antisense oligonucleotides into liver hepatocytes[J].Nucleic Acids Research, 2017(21):12388.DOI:10.1093/nar/gkx960.

