

Abstract


Targeted Protein Degradation (TPD) is considered one of the new drug discovery strategies that could address previously "undruggable" disease targets through the ubiquitin-proteasome system (UPS). Although the current Target Protein Centric approach has been proven successful in many cases, especially with well-validated E3 ligases such as CRBN and VHL, there is still a strong need to identify potential target protein degraders (PROTACs, molecular glues, etc.) through both target-centric and target-agnostic approaches. This would not only explore and identify potential tissue-specific novel E3 ligases but also greatly expand the range of potential degradable disease targets.

In this study, we first developed a molecular glue-focused library consisting of approximately 5,000 compounds through various medicinal chemistry evaluations. A GSPT1-HiBit cell line and a GSPT1-CRBN-KO counter cell line were developed for primary hit screening purposes. More than 100 potential hits were identified through the initial screening, which included both CRBN-dependent and CRBN-independent hits. A GSPT1-based proximity labeling cell line was then constructed for further target deconvolution purposes via proteomics.


Additionally, the synthetic-lethality mechanism can achieve specific killing of tumors with specific mutation background, laying the foundation for the development of precision treatment strategies based on molecular adhesives. We tried *in vitro* and *in vivo* synthetic-lethality screen to explore the sensitivity mechanism of liver cancer to PARP1 inhibition. We plan to future use the Molecular glue-focused library to identify novel degraders of PARP1 and other co-lethal genes, then combined both transcriptomics and proteomics assays to evaluate potential related E3 ligase. To sum up, our study demonstrates that integrated synthetic-lethality target discovery and target deconvolution strategies could be useful for TPD drug, realizing a closed loop from mechanism discovery to intervention strategies.

Synthetic Lethal Target Discovery Design


A. Introduction of synthetic-lethality strategy



B. Experiment design of synthetic-lethality screen




C. Workflow of synthetic-lethality target discovery




The **synthetic-lethality** mechanism can achieve specific killing of tumors with specific mutation background!
 A **whole-genome** library of mouse with **barcodes**² (~80k) for parallel *in vivo* and *in vitro* screen in liver cancer

Synthetic Lethal Target Discovery Process

A. Tumor volume cure for *in vivo* screen




B. Biomarker validation of screen samples



The changes of biomarker **PAR** indicated the PARP1 inhibit effect caused by Genotype-KO or Saruparib treatment.


Synthetic Lethal Target Discovery Result




The screening results identified **parp1** co-lethal genes that had been reported in other tumor models: **Dot1l** (Ovarian carcinoma³), **Apex1** (Gastric cancer⁴), **Pten** (Prostate cancer⁵, colorectal cancer⁶), **Tsc2** (Leiomyoma⁷), and also genes have not been reported before: **Map4k4**, **Ccm2**, **Zmynd8**...

Molecular Glue-focused Library Screen


A. Dose-response and uniformity test of ref compound




B. Z prime > 0.5 for all screening plates in the HTS campaign



C. Overview of tested compound




D. Dose-response of hit compounds across both target and counter cell lines




Nano-Glo® HiBIT assay at 384 well format was developed with robust performance.
 Screened **4.8K** compounds with hit rate **2.92%**.
 Hits identified over the counter screen include both **CRBN-dependent** and **CRBN-independent** degraders.

Biotin Ligase-based Proximity Labeling


A. Proximity Labeling E3 deconvolution strategy⁸



B. Biotinylation and pull down



C. LC-MS for E3 identification and ranking



CRBN could be captured and detected after GSPT1 degrader treatment for 4 hrs.
 The LC-MS analysis revealed multiple E3 ubiquitin ligases with a fold change exceeding 1.3.

Conclusion and Future Perspective

Conclusion

- Utilizing the construction of counter cell lines combined with high-throughput screening, CRBN-dependent and CRBN-independent degraders were successfully identified from the WuXi MGD focused library.
- Successfully established proximity labeling system combined with proteomics strategies to facilitate the deconvolution of E3 ligase.
- The iBar technology was utilized to achieve whole-genome *in vivo* screening, expanding the range of reported synthetic-lethality genes of PARP1 to liver cancer and identifying new synthetic-lethality genes.

Future Perspective

- We plan to future use the Molecular glue-focused library to identify novel degraders of PARP1 and other co-lethal genes, then combined both transcriptomics and proteomics assays to evaluate potential "degradable" targets.
- Potential candidate E3 ligases could be validated using a newly established array-based expression manipulation strategy.

References

- Linchong et al, 2021, Metabolic reprogramming and epigenetic modifications on the path to cancer, Protein Cell.
- Shiyu et al, 2019, Guide RNAs with embedded barcodes boost CRISPR-pooled screens.
- Liu, Chaohua et al, 2024. "PARP1-DOT1L transcription axis drives acquired resistance to PARP inhibitor in ovarian cancer." Molecular Cancer.
- Afzal, Hifza et al, 2019. "PARP1: A potential biomarker for gastric cancer." Pathology - Research and Practice.
- Bo et al, 2026, Suppression of PARP1 enhances PTEN mRNA therapy in castration-resistant prostate cancer by glycolysis disruption, Molecular Therapy Oncology.
- Mendes-Pereira et al, 2009. Synthetic lethal targeting of Pten mutant cells with parp inhibitors. Embo Molecular Medicine.
- Shikshya et al, 2022. ETV2 regulates PARP-1 binding protein to induce ER stress-mediated death in tuberlin-deficient cells. Life Science Alliance.
- Xiao et al, 2024, targeted protein degradation current and emerging approaches for E3 ligase deconvolution.