



Dianthus Platform | STING

Ready to Go Assays • Customizable • Cost-Effective • Fast Turnaround Time

Combining Ready-To-Go Assays with Dianthus TRIC Technology to Accelerate the Discovery of STING Modulators



Temperature-Related Intensity Change Technology

Sensing TRIC
Temperature Related Intensity Change



- Low sample consumption
- High throughput screening
- In solution binding affinities

2'3'-cGAMP binding to STING




Contact us

Technical Experts:

CreLux CR@wuxiapptec.com

DB_Service@wuxiapptec.com



wuxibiology.com



Search for more
protein info



Find us via
LinkedIn