


Bioconjugation Service


Bioconjugation in Early Drug Discovery

Target Identification


Biotin conjugation



PAL group conjugation




Fluorescence labeling




Hit Finding


Tb/Eu labeling for HTS



Enzymatic FRET substrate




Oligonucleotide coding




Efficiency Optimization


Antibody-based conjugation



Peptide-based conjugation



LNP-based conjugation



More than 20 years of organic synthesis experience from WuXi Chemistry.



Integrated screening platform, biochemical & biophysics modalities from WuXi Biology.



Experienced in providing diverse bioconjugation services throughout early discovery stage.

Conjugation Types

Macro-Molecules



Protein



Antibody



Peptide



Oligonucleotide




Beads

Linkers



Payloads




Toxin



Antigen



Si-RNA



Fluorophore

Chemistries & Linkers

Lysine

- Amide
- Carbamate
- Urea

Cysteine / Cystine

- Thioether
- Disulfide
- SMG


Non-natural Modification

- Click
- PAC
- Michael addition
- Hydrazone


Tags

- His-tag – NTA
- Halo-tag – Chloroalkane
- Snap-tag – OBG


Cleavable linkers




Hydrazone linker



Disulfide linker




Peptide linker




Glucuronide linker


Non-cleavable linkers




PEG linker




Thioether linker



Peptide linker




SGM linker



Triazole linker


Showcases and QC Methods

Antibody-Drug Conjugate (ADC) and DAR Measurement




ADC conjugation can be achieved by multiple biocompatible chemistries. DAR is determined by LCMS.

Cyclic Peptide-Oligomer-FAM Conjugate (cPOFC)




The conjugate is connected by click chemistry, QC by LCMS.




Cellular uptake of oligonucleotide is significantly increased by conjugating a cell penetration peptide.

Oligonucleotide-Drug Conjugate (ODC)



Cleavable Val-Cit-PAB linker can be quantitatively cleaved by cathepsin B and analyzed by LCMS.

Protein Oligonucleotide Conjugate (POC)



The conjugate is connected by NTA-His tag and analyzed by aSEC.



Contact Us

BD USA: mahnaz_arjomand@wuxiapptec.com
BD EU/UK/Israel: dave_madge@wuxiapptec.com

BD China: xu_longji@wuxiapptec.com
BD Japan: fumio_ito@wuxiapptec.com
BD Korea: sycho@wuxiapptec.com